



## BioFX™ Phosphate Buffered Saline/Casein Block and Diluent-5X Concentrate

Product insert

### Product number:

PBSC-0100-01 (100 mL)  
PBSC-1000-01 (1000 mL)  
Custom Sizes Available

### Product use:

BioFX™ Phosphate Buffered Saline/Casein Block and Diluent-5X Concentrate minimizes non-specific interactions of the secondary reagents with each other and the primary solid phase binding sites. The high blocking efficiency of casein is due to its content of small molecular weight proteins. The solution is ready for dilution to specific application requirements.

### Product stability, storage and specifications:

|                   |   |
|-------------------|---|
| Product stability | 2 years from the date of manufacture.   |
| Storage           | 2-8°C   |
| Specification     | This product should be clear and colorless in appearance. To minimize the risk of contamination, a fresh working solution should be prepared each day or store unused working solutions at the recommended storage temperature. |

### Recommendations for use:

- 1) Dilute 1 mL of PBS/Casein Block and Diluent-5X Concentrate in 4 mL reagent quality water to produce a 1X working solution.
- 2) Allow the 1X working solution to equilibrate to room temperature (25°C) prior to use.
- 3) To use as a blocking reagent, completely cover the target surface with the 1X working solution. To use as a diluent, dilute the antigen/antibody to a desired concentration in the 1X working solution and completely cover the target surface as needed.
  - An overnight incubation step at 2-8°C following the coating process will yield the best results.

For technical assistance, email [ivdtechsupport@surmodics.com](mailto:ivdtechsupport@surmodics.com)

### Related products:

|  |
|--|
| <b>In-Solution Protein Stabilizers &amp; Diluents:</b> |
| StabilZyme™ HRP Conjugate Stabilizer (SZ02)            |
| StabilZyme™ SELECT Stabilizer (SZ03)                   |
| StabilZyme™ NOBLE Stabilizer (SZ04)                    |
| StabilZyme™ Protein Free Stabilizer (SZPF)             |
| Surmodics™ Assay Diluent (Protein-Free) (SM01)         |
| MatrixGuard™ Diluent (SM02)                            |



**BioFX™ Phosphate Buffered Saline/Casein  
Block and Diluent-5X Concentrate**

Product insert

|   |
|---|
| <b>Blockers/Stabilizers:</b>  |
| StabilGuard™ Immunoassay Stabilizer-BSA-Free (SG01)                     |
| StabilCoat™ Immunoassay Stabilizer (SC01)                               |
| StabilBlock™ Immunoassay Stabilizer (ST01)                              |
| <b>Substrates:</b>  |
| BioFX™ TMB Super Sensitive One Component HRP Microwell Substrate (TMBS) |
| BioFX™ TMB One Component HRP Microwell Substrate (TMBW)                 |
| BioFX™ TMB Conductivity One Component HRP Microwell Substrate (TMBC)    |
| BioFX™ TMB Slow Kinetic One Component HRP Microwell Substrate (TMSK)    |
| BioFX™ TMB Extended Range One Component HRP Microwell Substrate (TMBX)  |
| BioFX™ TMB Enhanced HRP Membrane Substrate (ESPM)                       |
| BioFX™ TMB One Component HRP Membrane Substrate (TMBM)                  |
| BioFX™ Chemiluminescent Ultra Sensitive HRP Membrane Substrate (CHMM)   |

**SURMODICS, BIOFX, STABILZYME, MATRIXGUARD, STABILCOAT, STABILGUARD and STABILBLOCK are trademarks of Surmodics, Inc. and/or its affiliates.**

These products are for research and/or further manufacture and are not intended for use in humans or therapeutic purposes. Sales are without any seller's warranty or representation, expressed or implied, by usage or otherwise; no claims beyond replacement of unacceptable material or refund of purchase price shall be allowed. All claims must be made within 30 days following date of delivery.