

Product number:

TBSC-0100-01 (100 mL)
TBSC-1000-01 (1000 mL)

Intended use:

BioFX® Tris Buffered Saline/Casein Block and Diluent-10X Concentrate minimizes non-specific interactions of the secondary reactants with each other and with the primary solid phase binding sites. The high blocking efficiency of casein is due to its content of small molecular weight proteins. The solution is ready for dilution for specific application requirements.

Product stability, storage and specifications:

Product stability	2 years from date of manufacture.
Storage	2-8°C
Specification	The solution should be clear and colorless to light yellow in appearance. To minimize the risk of contamination, a fresh working solution should be prepared each day or store unused working solutions at the recommended storage temperature.

Recommendations for use:

- 1) Dilute 1 mL TBS Casein Block/Diluent-10X Concentrate in 9 mL reagent quality water to produce a 1X working solution.
- 2) Allow the 1X working solution to equilibrate to room temperature (25°C) prior to use.
- 3) To use as a blocking reagent, completely cover the target surface with the 1X working solution. To use as a diluent, dilute the antigen/antibody to a desired concentration in the 1X working solution and completely cover the target surface as needed.
 - An overnight incubation step at 2-8°C following the coating process will yield the best results.

Additional considerations:

TBS, pH 7.6, is a poor buffer below pH 7.5. The pH can fall more than 1 pH unit in warming from 0-37°C. Tris has the disadvantages of its reactivity as a primary amine and its appreciable solubility in organic solvents, which can lead to an accumulation in the biological phase of a reaction system. To avoid complex formation with ionic species such as calcium and magnesium in blood, tris buffers are preferable over phosphate buffers.

For technical assistance, email jvdttechsupport@surmodics.com

Related products:

BioFX® Phosphate Buffered Saline/Casein Block and Diluent – 5X Concentrate (PBSC)
BioFX® Chemiluminescent Ultra Sensitive HRP Membrane Substrate (CHMM)
BioFX® TMB Enhanced HRP Membrane Substrate (ESPM)

All products listed are for research use or for further manufacturing into *in vitro* diagnostic reagents. The products are not intended for use in humans or animals. Sales are without any seller's warranty or representation, expressed or implied, by usage or otherwise; no claims beyond replacement of product not meeting specifications or refund of purchase price shall be allowed. All claims must be made within 30 days following date of delivery.