

BBrMV
Banana bract mosaic virus
Catalog numbers: ECA 24600, SRA 24600

Technical Note:

There has been a change in the Agdia *Banana bract mosaic virus* (BBrMV) ELISA test. This change affects the enzyme conjugate diluents only. The BBrMV test now uses RUB6 buffer for the enzyme conjugate diluent instead of ECI. This diluent change will improve the OD values observed on negative tissue samples.

If you would like to use RUB6 with previously purchased reagents, Agdia will warranty the use of RUB6 with BBrMV reagents purchased from us on or after 8/29/2015.

If you have any questions about this product, please contact Agdia by phone (800-622-4342 or 574-264-2014) or email (info@agdia.com).