

# Tech Note: Protocol and Storage Change

ISK 49700 • InvictDetect™ ImmunoStrip®

---

## Technical Note

The extraction buffer for the InvictDetect ImmunoStrip has changed from AEB to AEB2. Ants must be homogenized in AEB2 buffer provided in kit.

For optimal results, each sample must contain five (5) ants. Testing more or fewer ants than recommended may result in invalid results.

The storage temperature for all kit components, including the ImmunoStrips, is now room temperature (18 - 30 °C).

For answers to your technical questions, please contact us at [techsupport@agdia.com](mailto:techsupport@agdia.com).

ImmunoStrip® is a registered trademark of Agdia, Inc.  
InvictDetect™ is a trademark of Agdia, Inc.



Agdia, Inc.  
52642 County Road 1  
Elkhart, IN 46514  
574-264-2014 / 800-622-4342  
[www.agdia.com](http://www.agdia.com) / [info@agdia.com](mailto:info@agdia.com)

n225  
Revised: 10/26/2020  
Page 1 of 1