

Tech Note: Test Improvement, Buffer Change

ISK/STX 49500 • *Tobacco etch virus* (TEV) ImmunoStrip®

Technical Note and Alert

Agdia is pleased to announce the availability of an improved ImmunoStrip for detection of Tobacco etch virus (TEV). The new test exhibits greater sensitivity and better flow characteristics than the previous version. Additionally, the new version uses SEB1 for the sample extraction buffer, making it easier to test the same sample extract for multiple pathogens.

Before using the test, please refer to the table below.

Item	Previous ImmunoStrip lot numbers	New ImmunoStrip lot numbers
ImmunoStrips	00004 and below	00005 and above
Sample Extraction Buffer	SEB9	SEB1
User Guide	m325.1	m175.5

For answers to your technical questions, please contact us at techsupport@agdia.com.

ImmunoStrip® is a registered Trademark of Agdia, Inc.



Agdia, Inc.
52642 County Road 1
Elkhart, IN 46514
574-264-2014 / 800-622-4342
www.agdia.com / info@agdia.com

n194
Revised: 05/01/2020
Page 1 of 1