

Bt-Cry34Ab1

Catalog numbers: 04500 and 11500

Technical Note:

Agdia is pleased to announce the availability of Bt-Cry34Ab1 ELISA test systems with a lower limit of detection for calculation of protein content.

Antibody components and plates from the new Agdia Bt-Cry34Ab1 tests are not interchangeable with antibody components and plates of the previous versions. Please refer to the table below for the listing of component compatibility.

Do not use previous test components with the new test components.

Please contact Agdia by phone (800-622-4342 or 574-264-2014) or email (info@agdia.com) if you have any questions about these test systems.

Item	Previous ELISA test 04500 lot numbers	New ELISA test 04500 lot numbers
Antibody-coated plates	00029 (solid)	00030 (solid)
Peroxidase enzyme conjugate	00309	00310
Instructions	m280.1	m280.2 and m280.3

Item	Previous ELISA test 11500 lot numbers	New ELISA test 11500 lot numbers
Antibody-coated plates	00026 (1X8)	00031 (1X8)
Peroxidase enzyme conjugate	00304	00305
Instructions	m281.2	m281.3