

# Content List: PathoScreen® Kits

## Complete Kit • Components and Sizes

### 96 Testwells

Component	Size/Number
Antibody-coated 96-well microtiter plates	1
Alkaline phosphatase enzyme conjugate	0.150 mL
RUB6 enzyme conjugate diluent, 1X ready to use	11 mL
PNP substrate tablets, 5 mg each	12 count
PNP substrate buffer, 5X concentrate	12 mL
Positive control	1
PBST buffer, 20X concentrate, 50 mL	3
TWEEN® 20	15 mL
General extract buffer 3, GEB3	24 g

### 288 Testwells

Component	Size/Number
Antibody-coated 96-well microtiter plates	3
Alkaline phosphatase enzyme conjugate	0.275 mL
RUB6 enzyme conjugate diluent, 1X ready to use	33 mL
PNP substrate tablets, 5 mg each	12 count
PNP substrate buffer, 5X concentrate	12 mL
Positive control	1
PBST buffer, 20X concentrate, 50 mL	5
TWEEN® 20	30 mL
General extract buffer 3, GEB3	48 g

### 480 Testwells

Component	Size/Number
Antibody-coated 96-well microtiter plates	5
Alkaline phosphatase enzyme conjugate	0.525 mL
RUB6 enzyme conjugate diluent, 1X ready to use	55 mL
PNP substrate tablets, 5 mg each	25 count
PNP substrate buffer, 5X concentrate	25 mL
Positive control	1
PBST buffer, 20X concentrate, 50 mL	7
TWEEN® 20	30 mL
General extract buffer 3, GEB3	48 g

For answers to your technical questions, please contact us at [techsupport@agdia.com](mailto:techsupport@agdia.com).

PathoScreen® is a registered trademark of Agdia, Inc.  
TWEEN® is a registered trademark of Croda International PLC.



Agdia, Inc.  
52642 County Road 1  
Elkhart, IN 46514  
574-264-2014 / 800-622-4342  
[www.agdia.com](http://www.agdia.com) / [info@agdia.com](mailto:info@agdia.com)