

Content List: PathoScreen® Kits

Complete Kit • Components and Sizes

96 Testwells

Component	Size/Number
Antibody-coated 96-well microtiter plates	1
Alkaline phosphatase enzyme conjugate	0.150 mL
ECl buffer, 5X concentrate	3 mL
PNP substrate tablets, 5 mg each	12 count
PNP substrate buffer, 5X concentrate	12 mL
Positive control	1
PBST buffer, 20X concentrate, 50 mL	3
TWEEN® 20	15 mL
General extract buffer, GEB	16.5 g

288 Testwells

Component	Size/Number
Antibody-coated 96-well microtiter plates	3
Alkaline phosphatase enzyme conjugate	0.275 mL
ECl buffer, 5X concentrate	6.5 mL
PNP substrate tablets, 5 mg each	12 count
PNP substrate buffer, 5X concentrate	12 mL
Positive control	1
PBST buffer, 20X concentrate, 50 mL	5
TWEEN® 20	30 mL
General extract buffer, GEB	33 g

480 Testwells

Component	Size/Number
Antibody-coated 96-well microtiter plates	5
Alkaline phosphatase enzyme conjugate	0.525 mL
ECl buffer, 5X concentrate	11 mL
PNP substrate tablets, 5 mg each	25 count
PNP substrate buffer, 5X concentrate	25 mL
Positive control	1
PBST buffer, 20X concentrate, 50 mL	7
TWEEN® 20	30 mL
General extract buffer, GEB	33 g

For answers to your technical questions, please contact us at techsupport@agdia.com.

PathoScreen® is a registered trademark of Agdia, Inc.
TWEEN® is a registered trademark of Croda International PLC.



Agdia, Inc.
52642 County Road 1
Elkhart, IN 46514
574-264-2014 / 800-622-4342
www.agdia.com / info@agdia.com