

DMO-RR-2A-1Ac ImmunoStrip®

Catalog number: 42900

Technical Note:

The buffer and user guide for the DMO-RR-2A-1Ac ImmunoStrip have changed.

Before using this test, please refer to the table below.

Item	Previous ImmunoStrip lot numbers	New ImmunoStrip lot numbers
ImmunoStrips	00019 and below	00020 and above
Sample extraction buffer	SEB4	PBST
Tissue	Leaf and Seed	Seed
User guide	m381	m381.1

Technical Support:

Please contact Agdia by phone (800-622-4342 or 574-264-2014) or email (info@agdia.com) if you have any questions about this test.