



BioFX™ 650 nm Liquid Stop Solution for TMB Microwell Substrates

Product insert

Product number:

LBSP-0100-01 (100 mL)
LBSP-1000-01 (1000 mL)
Custom Sizes Available

Product use:

BioFX™ 650 nm Liquid Stop Solution for TMB Microwell Substrates is a propriety formulation used to stop TMB substrate reactions in the blue spectral range.

Product stability, storage and specifications:

Product stability	3 years from the date of manufacture.
Storage	Storage recommendations are 2-25°C.
Specification	Blue color is read at 620-650 nm.

Recommendations for use:

Equal volumes of substrate and stop solution should be used to stop the TMB reaction.

For technical assistance, email ivdtechsupport@surmodics.com

Related products:

Substrates:
BioFX™ TMB Super Sensitive One Component HRP Microwell Substrate (TMBS)
BioFX™ TMB One Component HRP Microwell Substrate (TMBW)
BioFX™ TMB Conductivity One Component HRP Microwell Substrate (TMBC)
BioFX™ TMB Slow Kinetic One Component HRP Microwell Substrate (TMSK)
BioFX™ TMB Extended Range One Component HRP Microwell Substrate (TMBX)
TMB Stop Solutions:
BioFX™ 450 nm Liquid Nova-Stop Solution for TMB Microwell Substrates (NSTP)
BioFX™ 450 nm Liquid Stop Solution for TMB Microwell Substrates (LSTP)
BioFX™ 450 nm Stop Reagent for TMB Microwell Substrates (STPR)
BioFX™ 650 nm Stop Reagent for TMB Microwell Substrates (BSTP)

SURMODICS and BIOFX are trademarks of Surmodics, Inc. and/or its affiliates.

These products are for research and/or further manufacture and are not intended for use in humans or therapeutic purposes. Sales are without any seller's warranty or representation, expressed or implied, by usage or otherwise; no claims beyond replacement of unacceptable material or refund of purchase price shall be allowed. All claims must be made within 30 days following date of delivery.