

Tech Note: Coating Antibody Information

SRA/PSA 44800 • *Cucumber mosaic virus*-serogroup II (CMV II) ELISA

Technical Note

Pathoscreen® kits for CMV II use plates coated with antibody for *Cucumber mosaic virus* (CMV). These plates capture all strains of CMV. The enzyme conjugate in the kit is specific for CMV II.

Reagent sets for CMV II include coating antibody for *Cucumber mosaic virus* (CMV). This coating antibody captures all strains of CMV. The enzyme conjugate in the reagent set is specific for CMV II.

For answers to your technical questions, please contact us at techsupport@agdia.com.

PathoScreen® is a registered Trademark of Agdia, Inc



Agdia, Inc.
52642 County Road 1
Elkhart, IN 46514
574-264-2014 / 800-622-4342
www.agdia.com / info@agdia.com

n29
Revised: 09/04/2019
Page 1 of 1